



Roche © 2009 All rights reserved		CPDG A/C COMBO EXPERT DM			
Revisions			For revision approval date See CN		
Rev	Chg #	Description	Date	Reviser	Approver
B	RDODC1357	Updated to support DM LCM Project	7/14/2009	Rick Wilson	James Parker
C	RDODC2675	Updated to support DM LCM2 Project	4/27/2010	Rick Wilson	James Parker
D	RDODC3690	Updated to support DM Brazilian Portuguese	1/21/2011	Bob Sabo	Chuck Bolam
E	RDODC6107	Updated to reflect updated strategy regarding third party commands.	9/21/2012	Bob Sabo	Chuck Bolam

ACCU-CHEK Combo Aviva, Performa, ACCU-CHEK Expert Aviva, Performa

Communications Protocol Developer's Guideline

**This document is provided as part of a license agreement only.
Use of the contents of this document without a license
agreement with Roche Diagnostics Operations is prohibited.**

Roche Diagnostics Operations makes no representations or warranties with respect to the contents of this documentation and specifically disclaims any implied warranties, including the implied warranties of merchantability and fitness for a particular purpose. In no case shall Roche Diagnostics Operations be liable for incidental or consequential damages.

Document is subject to change without notification.

	Originators	Date	Release Doc. No. RDODC629
Originator	Rick Wilson	03/17/2009	
Engineer	Ray Strickland	03/17/2009	
Engineer			
Approved	Karl Werner	See CN	Ref Roche Part No. --

This document contains proprietary information and is loaned in confidence subject to return upon demand on the express condition that it will not be used in any way not authorized by Roche.

Type
SPEC

Document Number	Rev
7005079	E



Roche	CPDG A/C COMBO EXPERT DM
© 2009 All rights reserved	

Table of Contents

1	PURPOSE	3
2	DISCLAIMER	3
3	DEFINITIONS	3
4	REFERENCES	3
4.1	Manuals	3
4.2	Cables	4
5	DOCUMENT LEGEND	5
6	ACCU-CHEK METER PROTOCOL	6
6.1	Communication Requirements and Serial Port Settings	6
6.2	Protocol Command Structure	6
6.2.1	Data Blocks	7
6.2.2	Multiple Data Blocks	7
6.2.3	Checksum Calculator	7
6.2.4	Command interruption	7
6.3	Normal Sequence of Commands for Data Extraction	8
7	COMMANDS REFERENCE	10
7.1	Example of Commands with Parameters	10
8	ACCU-CHEK METER COMMANDS	11
8.1	Connect – <CAN>	11
8.2	Power Down – [1D]	11
8.3	Read & Clear Status – [0B]	11
8.4	Configuration – [43] or ‘C’	12
8.5	Instrument Name – [49] or ‘I’	15
8.6	Read Setup – [53] or ‘S’	16
8.7	Change Setup – [0C]	37
8.8	Enable/Disable Timers – [5A] or ‘Z’	42
8.9	Obtain Number of Results – [60]	44
8.10	Send Results from Start to End – [61] or ‘a’	44
8.11	Reset Results Memory – [52] or ‘R’	48
9	STATUS REGISTER VALUES	49

Roche Confidential Document

Released

**Company
Confidential**

Page 2 of 49

Document Number	Rev
7005079	E



Roche © 2009 All rights reserved	CPDG A/C COMBO EXPERT DM
--	---------------------------------

1 Purpose

The purpose of this document is to provide the minimum communication protocol for external developers to extract results from an ACCU-CHEK® blood glucose meter. The audiences of this document are software development professionals who are interested in developing an interface between Roche Diagnostics ACCU-CHEK meter and an external host.

This document only provides information for Roche Diagnostics ACCU-CHEK® Combo and ACCU-CHEK® Expert meters. Other Roche Diagnostic meters are addressed in separate CPDG documents.

2 Disclaimer

The developer using this document to extract data from a Roche Diagnostics meter assumes all responsibility for conforming to this standard. Data extracted from a Roche Diagnostics meter is not to be manipulated by an external process. Data should be used only for transmission to a central repository or data analysis system. Roche Diagnostics is not responsible for inaccurate or misdiagnosis of a patient condition due to mishandling of data extracted from a Roche Diagnostics blood glucose meter.

3 Definitions

bG – Abbreviation for **blood Glucose**.

cG – Term used to refer to test run with control solution.

Host – Device communicating to meter, such as a PC.

IR – Infrared light

Monitor – A blood Glucose monitor, referred to as a meter.

Number – Any sequence of numeric characters ('0'-'9').

Date – A string of fixed length ASCII characters that form a date using YYMMDD format.

Time – A string of fixed length ASCII characters that represent a time using HHMMSS format.

Function – A single ASCII character used as a sub-division of a command.

Fixed Hex Word – Any string of 4 ASCII characters that a hexadecimal number can be formed.

Fixed Hex Double Word – Any string of 8 ASCII characters that a hexadecimal number can be formed.

Hex Word – Any string of 1 to 4 ASCII characters that a hexadecimal number can be formed.

Fixed Hex Byte – Any string of 2 ASCII characters that a hexadecimal number can be formed.

Hex Byte – Any string of 1 to 2 ASCII characters that a hexadecimal number can be formed.

String – Any string of printable ASCII characters.

Float – Any String of 1 to 4 characters that a floating point number can be formed from. A decimal point is included in the string.

Signed Number – Any sequence of numeric characters ('0'-'9') preceded by a "+" or a "-".

4 References

4.1 Manuals

ACCU-CHEK Aviva Combo User Manual
ACCU-CHEK Performa Combo User Manual
ACCU-CHEK Aviva Expert User Manual
ACCU-CHEK Performa Expert User Manual

**Company
Confidential**

Page 3 of 49

Document Number	Rev
7005079	E



Roche © 2009 All rights reserved	CPDG A/C COMBO EXPERT DM
--	---------------------------------

4.2 Cables

ACCU-CHEK Infrared Cable (IR-210B Infrared Adapter by TekRam Technology)
Roche Diagnostics Cat. No. 3183408

ACCU-CHEK Universal Cable: Roche Diagnostics Cat. No. 03062678001
This document does not provide information for connection to the Roche Diagnostics USB interface cable.

Roche Confidential Document
Released

Company Confidential

Document Number 7005079	Rev E
-----------------------------------	-----------------



Roche © 2009 All rights reserved	CPDG A/C COMBO EXPERT DM
--	---------------------------------

5 Document Legend

Legend:

{ } Hex value represented using ASCII
< > Serial Control Character
[] Hex value
→ From host to instrument
← From instrument to host
^ Control character – Press character on keyboard while holding the CTRL key down.

Example

{4A}
<CR>
[4A]

Description

Two bytes, '4' and 'A', are placed into the data stream.
Carriage return is placed into stream.
4A hex gets placed into stream.

Serial Control Characters:

<STX> = [02]	<ACK> = [06]	<LF> = [0A]	<ETX> = [03]
<TAB> = [09]	<NAK> = [15]	<CR> = [0D]	<EOT> = [04]
<CAN> = [18]			

Key Points

Input and Output parameters are <TAB> delimited.

All BG results are transmitted in mg/dL units. If mmol/L units are desired, then they **must** be converted using the formula: 1 mmol/L = 18.02 mg/dl.

Roche Confidential Document

Released

**Company
Confidential**

Page 5 of 49

Document Number
7005079

Rev
E



Roche © 2009 All rights reserved	CPDG A/C COMBO EXPERT DM
--	---------------------------------

6 ACCU-CHEK Meter Protocol

6.1 Communication Requirements and Serial Port Settings

The ACCU-CHEK meter serial communications are at **9600** baud with **8** data bits, **1** start bit, **1** stop bit and **NO** parity checking. The ACCU-CHEK meter communicates via Infrared (IR). In order for a Host to communicate with the meter, an IR Cable must be attached to the PC. Cables for this function are available from Roche Diagnostics. The ACCU-CHEK meter **must be in Communication mode** to allow serial communications. From the "My Data" menu screen select "Data Transfer" to get the meter into Communication mode. The automatic power off time is 2 minutes.

The protocol shall only support half-duplex communications. When switching from receiving to transmitting, the meter shall insert a 10ms delay before transmitting. When switching from transmitting to receiving, the meter shall expect 10ms delay from the host.

When the meter is receiving data frames, the meter requires a two-millisecond delay between each character. When the meter is transmitting data frames, an inter-character delay may be present.

6.2 Protocol Command Structure

Each command consists of at least a command (one character) and a command terminator. The number of parameters for each command and the allowed parameter range are specified in the command reference.

The meter will echo back each command character and each single parameter character to increase the reliability of the communication protocol as soon as it is received. The command terminator is not echoed back.

Commands can be accepted or rejected by the meter. If the meter accepts the command sent from the host, it will answer with an <ACK> directly after receiving the command terminator.

The meter shall reject commands by sending a <NAK> (after receiving the command terminator) if one of the following cases happened:

- the meter receives an unknown command
- the meter receives unexpected characters within the parameters
- the meter receives invalid number of parameters for the command
- one of the parameters is outside its allowed range
- the meter is in an error state, i.e., the status register is not zero

After receiving an accepted command and answering with <ACK> the firmware starts to execute the command.

After executing the command, the meter will send another <ACK or NAK> indicating successful/faulty completion of the command. Any problems while executing the command shall lead to storing the appropriate error number and answering with a <NAK>, otherwise the answer shall be an <ACK>.

If an error occurs in communication mode because of a wrong command or during the execution of the command, the firmware shall send a <NAK> as command response to the host and enter an error state. If the meter is in an error state, it will reject all commands by sending a <NAK> after receiving the command terminator until the error status is read out and cleared. The 'Read and Clear' command is the

**Company
Confidential**

Document Number	Rev
7005079	E



Roche © 2009 All rights reserved	CPDG A/C COMBO EXPERT DM
--	---------------------------------

only command accepted by the meter if it is in an error state. See Section 9 Status Register Values for definition of error codes.

6.2.1 DATA BLOCKS

Upload and download data shall be transferred in data blocks controlled by <STX> and <ETX or EOT> characters. These data blocks will have the following format:

<STX>{# of bytes}<TAB>Data<TAB>{CRC}<ETX or EOT>

STX: STX is the start of packet indicator.
of bytes: This is a two-byte ASCII string containing the number of data bytes, including the Tabs, as a hexadecimal number.
TAB: Used for Field delimitation.
Data: A stream of ASCII characters
TAB: Used for Field delimitation.
CRC: This is a two-byte ASCII string containing the 8-bit checksum value as a hexadecimal number.
ETX: ETX is used if more data blocks will be sent.
EOT: EOT is used if this is the last data block to be sent.

6.2.2 MULTIPLE DATA BLOCKS

In upload mode, the host has to send an <ACK or NAK> after each received data block. If the host sends, other characters, the meter aborts upload and sets the status to [FD], Aborted Command. In case of NAK the meter has to repeat the last data block. The host is notified that the meter has no more data blocks to send by the EOT terminator.

6.2.3 CHECKSUM CALCULATOR

To assure a reliable data transfer process, data blocks are protected by an 8-bit checksum (CRC). The checksum (CRC) is computed by bitwise XORing the data bytes with the previous checksum value. The initial value is [6E]. Only the <TAB>s and the data bytes are included in the checksum calculation, <STX> and # of bytes are not.

	STX	Length		TAB	Data					TAB	CRC		ETX/EOT
Data	<STX>	0	7	<TAB>	A	v	i	v	a	<TAB>	2	7	<EOT>
(Hex)	[02]	[30]	[37]	[09]	[41]	[76]	[69]	[76]	[61]	[09]	[32]	[37]	[04]
CRC	Initial value is [6E]			[67]	[26]	[50]	[39]	[4F]	[2E]	[27]			

6.2.4 COMMAND INTERRUPTION

The host as the communication master has the possibility to interrupt the meter receiving commands by sending a <CAN> character. This cancel command might be necessary e.g. after receiving an unexpected byte echo or wanting to abort a lengthy results download.

In this case the meter shall:

- Throw away all received command/parameter characters,

**Company
Confidential**

Document Number	Rev
7005079	E



Roche © 2009 All rights reserved	CPDG A/C COMBO EXPERT DM
--	---------------------------------

- Answer a <NAK> to the host,
- Wait for new commands from the host.

At this time, the Host will have to issue the Read and Clear Status Command before issuing any other commands.

The meter is capable of sending a <CAN>. If this happens, the meter will immediately send a <NAK> following the <CAN> and then return to command processing mode waiting for the next command.

6.3 Normal Sequence of Commands for Data Extraction

Example of initialization sequence required by ACCU-CHEK Infrared Cable*

Host Sends	Description
Set DTR := True Set RTS := False Wait (50 us) Send Control Byte = [14] Set RTS := True Set DTR := True; Wait (50 us);	[14] sets baud rate to 9600 and sets output pulse to 1.6 ms. See device technical spec for more information.

***No initialization sequence is needed for use with the ACCU-CHEK Universal Cable. However, it may be necessary to send a different initialization sequence if a different IR device (dongle) is used.**

Example of communication required for data extraction.

Host Sends	Meter Replies	Description
<CAN>	<NAK>	Initial communications
[0B] <CR> <ACK>	[0B] <ACK> <STX>Response<EOT> <ACK>	Read and clear Status
[43] [09] [31] <CR> <ACK>	[43] [09] [31] <ACK> <STX>Response<EOT> <ACK>	Read Software Version command
[43] [09] [34] <CR> <ACK>	[43] [09] [34] <ACK> <STX>Response<EOT> <ACK>	Read Model number command
[43] [09] [33] <CR> <ACK>	[43] [09] [33] <ACK> <STX>Response<EOT> <ACK>	Read Serial number command
[53] [09]	[53] [09]	Read bG units command

**Company
Confidential**

Document Number
7005079

Rev
E



Roche © 2009 All rights reserved	CPDG A/C COMBO EXPERT DM
--	---------------------------------

[33] <CR> <ACK>	[33] <ACK> <STX>Response<EOT> <ACK>	
[60] <CR> <ACK>	[60] <ACK> <STX>Response<EOT> <ACK>	Get number of results stored
[61] [09] [xx] [09] [yy] <CR> <ACK>	[61] [09] [xx] [09] [yy] <ACK> <STX>Response<ETX> or <EOT> <ACK>	Extract meter records from memory
[1D] <CR>	[1D] <ACK> <ACK>	Power down command

* Response includes # of bytes, <TAB>, data, <TAB> and CRC

* xx is the first requested result record and yy is the last requested result record.

Roche Confidential Document
Released

Company Confidential

Document Number 7005079	Rev E
-----------------------------------	-----------------



Roche © 2009 All rights reserved	CPDG A/C COMBO EXPERT DM
--	---------------------------------

7 Commands Reference

7.1 Example of Commands with Parameters

Example of Commands with parameters			
Command Character This section contains the command character.	<i>This section contains a general description of the command.</i>		
This section contains the command character plus, if present, the command sub-function. Both the command and the function are one character separated by a <TAB>.	Input/Output Name This column contains the name of the given parameter.	Input/Output Type This column contains the data type information. The data types are defined in Section 3 Definitions.	Input/Output Range This column contains limits or range of the given parameter.
	This section typically contains an example of the command, plus additional information such as notes and sources.		

Roche Confidential Document
Released

Company Confidential

Document Number 7005079	Rev E
-----------------------------------	-----------------



Roche © 2009 All rights reserved	CPDG A/C COMBO EXPERT DM
--	---------------------------------

8 ACCU-CHEK Meter Commands

8.1 Connect – <CAN>	
Command Character <CAN>	The <CAN> command is used to initialize communication with the meter. <CAN> may also be used to interrupt commands in process.
	Example: Turn on the meter. →<CAN> ←<NAK> Note: The host must issue the Read and Clear Status command ([0B]) to clear the NAK response. If the <CAN> command is issued again, a <NAK> character is returned.

8.2 Power Down – [1D]	
Command Character [1D]	The [1D] command shall power down the meter to “off” state.
	Example: Turn off the meter, power down in “off” state. →[1D] ←[1D] →<CR> ←<ACK> ←<ACK>

8.3 Read & Clear Status – [0B]			
Command Character [0B]	The [0B] command shall read and clear the meter status.		
	Output Name Status Register	Output Type Fixed Hex Word	Output Range 0000H – FFFFH
	Example: Initial instrument communications established →[0B] ←[0B] →<CR> ←<ACK> ←<STX>{LEN}<TAB>{ASCII Data}<TAB>{CRC}<EOT> →<ACK> ←<ACK> Note: The status register values are defined in the “Status Register Values” section of this document. See section 9 for error code descriptions.		

**Company
Confidential**

Page 11 of 49

Document Number
7005079

Rev
E



Roche © 2009 All rights reserved	CPDG A/C COMBO EXPERT DM
--	---------------------------------

8.4 Configuration – [43] or ‘C’			
Command Character [43] or ‘C’	The [43] or ‘C’ commands shall be used to read information about the meter’s configuration.		
C<TAB>1 — Read software version number	The C<TAB>1 command shall return the meter software version number.		
	Output Name Software Version	Output Type String	Output Range 1 – 7 characters
	Example: Version number = 4.11 →C ←C →<TAB> ←<TAB> →1 ←1 →<CR> ←<ACK> ←<STX>{LEN}<TAB>4.11<TAB>{CRC}<EOT> →<ACK> ←<ACK>		
C<TAB>2 — Read meter hardware version number	The C<TAB>2 command shall return the meter’s hardware version.		
	Output Name Hardware Version	Output Type String	Output Range 1 – 8 characters
	Example: Version number = DM.01 →C ←C →<TAB> ←<TAB> →2 ←2 →<CR> ←<ACK> ←<STX>{LEN}<TAB>DM.01<TAB>{CRC}<EOT> →<ACK> ←<ACK>		

Company Confidential

Document Number 7005079	Rev E
-----------------------------------	-----------------

Roche Confidential Document
Released



Roche © 2009 All rights reserved	CPDG A/C COMBO EXPERT DM
--	---------------------------------

8.4 Configuration – [43] or ‘C’

C<TAB>3 — Read instrument serial number	The C<TAB>3 command shall return the meter's serial number.		
	Output Name Serial Number	Output Type String	Output Range 1 – 11 characters
	Example: Serial number = 70000001 →C ←C →<TAB> ←<TAB> →3 ←3 →<CR> ←<ACK> ←<STX>{LEN}<TAB>70000001<TAB>{CRC}<EOT> →<ACK> ←<ACK>		
C<TAB>4 — Read model number	The C<TAB>4 command shall return the meter's model number.		
	Output Name Model Number	Output Type String	Output Range 1 – 7 characters
	Example: Model number = 535 →C ←C →<TAB> ←<TAB> →4 ←4 →<CR> ←<ACK> ←<STX>{LEN}<TAB>535<TAB>{CRC}<EOT> →<ACK> ←<ACK>		

**Company
Confidential**

Document Number	Rev
7005079	E



Roche © 2009 All rights reserved	CPDG A/C COMBO EXPERT DM
--	---------------------------------

8.4 Configuration – [43] or ‘C’

C<TAB>5 — Read transmit time (in sec) for max. records @ 9600 baud	The C<TAB>5 command shall return the time to transmit the maximum number of records at 9600 baud.		
	Output Name Max. Transmit Time	Output Type Number	Output Range 0 – 65535
	Example: Transmit time = 91 seconds →C ←C →<TAB> ←<TAB> →5 ←5 →<CR> ←<ACK> ←<STX>{LEN}<TAB>091<TAB>{CRC}<EOT> →<ACK> ←<ACK> Note: The transmit time reported is based on a theoretical model and will vary from the actual time.		
C<TAB>6 — Read maximum number of records that can be stored	The C<TAB>6 command shall return the sum of the maximum number of control records and bG records that can be stored.		
	Output Name Max. Number of Records	Output Type Number	Output Range 0 – 1100
	Example: Maximum number of records = 1100 →C ←C →<TAB> ←<TAB> →6 ←6 →<CR> ←<ACK> ←<STX>{LEN}<TAB>1100<TAB>{CRC}<EOT> →<ACK> ←<ACK>		

**Company
Confidential**

Document Number	Rev
7005079	E

Roche Confidential Document
Released



Roche © 2009 All rights reserved	CPDG A/C COMBO EXPERT DM
--	---------------------------------

8.4 Configuration – [43] or 'C'

C<TAB>b — Determine Type of meter (Combo or Expert).	The C<TAB>b command shall return the 0 if meter is an ACCU-CHEK® Expert Meter and return 1 if the meter is an ACCU-CHEK® Combo meter.		
	Output Name Meter Type	Output Type Number	Output Range 0 – 1 0 – Expert 1 – Combo
	Example: Meter Type = Expert →C ←C →<TAB> ←<TAB> →b ←b →<CR> ←<ACK> ←<STX>{LEN}<TAB>0<TAB>{CRC}<EOT> →<ACK> ←<ACK>		

8.5 Instrument Name – [49] or 'I'

Command Character [49] or 'I'	The [49] or 'I' command shall be used to read the instrument name.		
	Output Name Name	Output Type String	Output Range 1 – 32 characters
	Example: Name = DM Aviva →I<CR> ←I<ACK> ←<STX>{LEN}<TAB>DM Aviva<TAB>{CRC}<EOT> →<ACK> ←<ACK> Note: The 'DM Aviva' in this example may vary by model number. This is an example only.		

**Company
Confidential**

Document Number	Rev
7005079	E



Roche © 2009 All rights reserved	CPDG A/C COMBO EXPERT DM
--	---------------------------------

8.6 Read Setup – [53] or ‘S’

Command Character [53] or ‘S’	The [53] or ‘S’ commands shall allow the host to extract the meter settings from the meter.		
S<TAB>1 — Read date	The S<TAB>1 command shall read the date.		
	Output Name Date	Output Type Date	Output Range 090101 – 311231 Date format is YYMMDD
	Example: Date = February 3, 2009 →S ←S →<TAB> ←<TAB> →1 ←1 →<CR> ←<ACK> ←<STX>{LEN}<TAB>090203<TAB>{CRC}<EOT> →<ACK> ←<ACK>		
S<TAB>2 — Read time	The S<TAB>2 command shall read the time.		
	Output Name Time	Output Type Time	Output Range 000000 – 235959 Time format is 24-hour (HHMMSS)
	Example: Time = 3:45 PM →S ←S →<TAB> ←<TAB> →2 ←2 →<CR> ←<ACK> ←<STX>{LEN}<TAB>154500<TAB>{CRC}<EOT> →<ACK> ←<ACK>		

**Company
Confidential**

Document Number
7005079

Rev
E

Roche Confidential Document
Released



Roche © 2009 All rights reserved	CPDG A/C COMBO EXPERT DM
--	---------------------------------

8.6 Read Setup – [53] or ‘S’

S<TAB>3 — Read units	The S<TAB>3 command shall read the display units.		
	Output Name Units	Output Type String	Output Range “mg/dl” or “mmol/l”
	Example: Units = mg/dl →S ←S →<TAB> ←<TAB> →3 ←3 →<CR> ←<ACK> ←<STX>{LEN}<TAB>mg/dl<TAB>{CRC}<EOT> →<ACK> ←<ACK>		
S<TAB>6 — Read Patient Hypo. Threshold	The S<TAB>6 command shall return the Hypo level.		
	Output Name Hypo. Threshold	Output Type Number (mg/dl)	Output Range 50 – 90
	Example: Hypo. Threshold = 66 mg/dl →S ←S →<TAB> ←<TAB> →6 ←6 →<CR> ←<ACK> ←<STX>{LEN}<TAB>66<TAB>{CRC}<EOT> →<ACK> ←<ACK> Note: The Output Range is 54 – 90 mg/dl for a meter that was configured as a mmol/l meter.		

**Company
Confidential**

Document Number	Rev
7005079	E



Roche © 2009 All rights reserved	CPDG A/C COMBO EXPERT DM
--	---------------------------------

8.6 Read Setup – [53] or ‘S’

S<TAB>7 — Read Selected Language	The S<TAB>7 command shall read the selected language.		
	Output Name Selected Language	Output Type Number	Output Range 1 – 7, A, C, D, F, G, H, I, J, M, N, O, Q 1 = German 2 = English 3 = French 4 = Spanish 5 = Italian 6 = Dutch 7 = Swedish A = Portuguese C = Norwegian D = Finnish F = Danish G = Czech H = Hungarian I = Slovenian J = Russian M = Polish N = Slovakian O = Brazilian Portuguese Q = Romanian
	Example: Language = English →S ←S →<TAB> ←<TAB> →7 ←7 →<CR> ←<ACK> ←<STX>{LEN}<TAB>2<TAB>{CRC}<EOT> →<ACK> ←<ACK>		

**Company
Confidential**

Document Number	Rev
7005079	E



Roche © 2009 All rights reserved	CPDG A/C COMBO EXPERT DM
--	---------------------------------

8.6 Read Setup – [53] or ‘S’

S<TAB>8 — Read date format	The S<TAB>8 command shall read the date format.		
	Output Name Date Format	Output Type Number	Output Range 3 3 = DDMMYY
	Example: Date Format = DDMMYY →S ←S →<TAB> ←<TAB> →8 ←8 →<CR> ←<ACK> ←<STX>{LEN}<TAB>3<TAB>{CRC}<EOT> →<ACK> ←<ACK> Note: The meter always return 3 for this command because only one date format is supported. (DDMMYY; for example: 01Jan09.)		
S<TAB>9 — Read time format	The S<TAB>9 command shall read the time format.		
	Output Name Time Format	Output Type Number	Output Range 1 – 2 1 = 24 2 = 12
	Example: Time Format = 12-Hour →S ←S →<TAB> ←<TAB> →9 ←9 →<CR> ←<ACK> ←<STX>{LEN}<TAB>2<TAB>{CRC}<EOT> →<ACK> ←<ACK>		

**Company
Confidential**

Document Number	Rev
7005079	E

Roche Confidential Document
Released



Roche © 2009 All rights reserved	CPDG A/C COMBO EXPERT DM
--	---------------------------------

8.6 Read Setup – [53] or ‘S’

S<TAB>H — Read Insulin Increment	The S<TAB>H command shall read the insulin increment setting.		
	Output Name Insulin Increment	Output Type Number	Output Range 0 – 2 0 = 0.1 1 = 0.5 2 = 1.0
	Example: Insulin Increment = 1.0 →S ←S →<TAB> ←<TAB> →H ←H →<CR> ←<ACK> ←<STX>{LEN}<TAB>2<TAB>{CRC}<EOT> →<ACK> ←<ACK>		

Roche Confidential Document
Released

Company Confidential

Document Number 7005079	Rev E
-----------------------------------	-----------------



Roche © 2009 All rights reserved	CPDG A/C COMBO EXPERT DM
--	---------------------------------

S<TAB>Q — Read Stored Languages	The S<TAB>Q command shall return the language IDs for the stored languages. The meter can contain from 1 to 4 languages in its non-volatile memory. This command returns the language ID for each language that is stored.		
	Output Name Language IDs	Output Type Number	Output Range 1 – 7, A, C, D, F, G, H, I, J, M, N, O, Q 1 = German 2 = English 3 = French 4 = Spanish 5 = Italian 6 = Dutch 7 = Swedish A = Portuguese C = Norwegian D = Finnish F = Danish G = Czech H = Hungarian I = Slovenian J = Russian M = Polish N = Slovakian O = Brazilian Portuguese Q = Romanian
	Example: The meter has the following languages stored: English, German, French, and no language stored in its 4th Language Section. →S ←S →<TAB> ←<TAB> →Q ←Q →<CR> ←<ACK> ←<STX>{LEN}<TAB>2<TAB>1<TAB>3<TAB><TAB>{CRC}<EOT> →<ACK> ←<ACK>		

**Company
Confidential**

Document Number 7005079	Rev E
-----------------------------------	-----------------



Roche © 2009 All rights reserved	CPDG A/C COMBO EXPERT DM
--	---------------------------------

8.6 Read Setup – [53] or ‘S’

S<TAB>T — Read Carb Units	The S<TAB>T command shall read the Carb Units setting.		
	Output Name Carb Units	Output Type Number	Output Range 0, 1, 4, 5 0 = Grams 1 = BE 4 = CC 5 = KE
	Example: Carb Units = Grams →S ←S →<TAB> ←<TAB> →T ←T →<CR> ←<ACK> ←<STX>{LEN}<TAB>0<TAB>{CRC}<EOT> →<ACK> ←<ACK>		
S<TAB>[6D] — Read Bolus Advice Status	The S<TAB>[6D] command shall return the bolus advice status.		
	Output Name Advice Status	Output Type Number	Output Range 0 – 1 (0 = disabled)
	Example: Advice Status = Disabled →S ←S →<TAB> ←<TAB> →[6D] ←[6D] →<CR> ←<ACK> ←<STX>{LEN}<TAB>0<TAB>{CRC}<EOT> →<ACK> ←<ACK>		

**Company
Confidential**

Page 22 of 49

Document Number
7005079

Rev
E



Roche © 2009 All rights reserved	CPDG A/C COMBO EXPERT DM
--	---------------------------------

S<TAB>[6E] — Read Key Lock Status	The S<TAB>[6E] command shall return the key lock status.		
	Output Name Key Lock Status	Output Type Number	Output Range 0 – 1 (0 = disabled)
	Example: Key Lock Status = Disabled →S ←S →<TAB> ←<TAB> →[6E] ←[6E] →<CR> ←<ACK> ←<STX>{LEN}<TAB>0<TAB>{CRC}<EOT> →<ACK> ←<ACK>		
S<TAB>[6F] — Read Setup Wizard Done Status	The S<TAB>[6F] command shall return the Setup Wizard Done status.		
	Output Name Setup Wizard Done Status	Output Type Number	Output Range 0 – 1 (0 = not done)
	Example: Setup Wizard Not Done →S ←S →<TAB> ←<TAB> →[6F] ←[6F] →<CR> ←<ACK> ←<STX>{LEN}<TAB>0<TAB>{CRC}<EOT> →<ACK> ←<ACK>		

**Company
Confidential**

Document Number 7005079	Rev E
-----------------------------------	-----------------

Roche Confidential Document
Released



Roche © 2009 All rights reserved	CPDG A/C COMBO EXPERT DM
--	---------------------------------

8.6 Read Setup – [53] or ‘S’

S<TAB>[70] — Read Bluetooth Control Status	The S<TAB>[70] command shall read the Bluetooth control status.		
	Output Name Bluetooth Control Status	Output Type Number	Output Range 0 – 1 (0 = Bluetooth off)
	Example: Bluetooth Control Status = Off →S ←S →<TAB> ←<TAB> →[70] ←[70] →<CR> ←<ACK> ←<STX>{LEN}<TAB>0<TAB>{CRC}<EOT> →<ACK> ←<ACK> Note: This command is for ACCU-CHEK® Combo meters only. ACCU-CHEK® Expert meters support this command and always return zero.		
S<TAB>[71] — Read Key Click Status	The S<TAB>[71] command shall return the key click status.		
	Output Name Key click Status	Output Type Number	Output Range 0 – 1 (0 = Key click disabled)
	Example: Key click Status = Disabled →S ←S →<TAB> ←<TAB> →[71] ←[71] →<CR> ←<ACK> ←<STX>{LEN}<TAB>0<TAB>{CRC}<EOT> →<ACK> ←<ACK>		

**Company
Confidential**

Page 24 of 49

Document Number
7005079

Rev
E



Roche © 2009 All rights reserved	CPDG A/C COMBO EXPERT DM
--	---------------------------------

S<TAB>[72] — Read Beeper Volume	The S<TAB>[72] command shall return the beeper volume setting.		
	Output Name Beeper Volume	Output Type Number	Output Range 0 – 3 0 = Off 1 = Low 2 = Medium 3 = High
	Example: Beeper Volume = Low →S ←S →<TAB> ←<TAB> →[72] ←[72] →<CR> ←<ACK> ←<STX>{LEN}<TAB>1<TAB>{CRC}<EOT> →<ACK> ←<ACK>		
S<TAB>[73] — Read Vibrate Status	The S<TAB>[73] command shall read the vibrate status.		
	Output Name Vibrate Status	Output Type Number	Output Range 0 – 1 (0 = Vibrate disabled)
	Example: Vibrate Status = Disabled →S ←S →<TAB> ←<TAB> →[73] ←[73] →<CR> ←<ACK> ←<STX>{LEN}<TAB>0<TAB>{CRC}<EOT> →<ACK> ←<ACK>		

Company Confidential

Document Number 7005079	Rev E
-----------------------------------	-----------------

Roche Confidential Document
Released



Roche © 2009 All rights reserved	CPDG A/C COMBO EXPERT DM
--	---------------------------------

S<TAB>[74] — Read bG Reminder	The S<TAB>[74] command shall return the bG Reminder Status, bG Reminder Threshold, and bG Reminder Duration for a selected bG Reminder Type.		
	Input Name bG Reminder Type	Input Type Number	Input Range 0 – 2 0 = After high bG 1 = After low bG 2 = After meal
	Output Name bG Reminder Status bG Reminder Threshold bG Reminder Time Duration	Output Type Number Number (mg/dl or grams) Time	Output Range 0 – 1 (0 = disabled) <i>Thres Min – Thres Max</i> For After high bG, the populated range is 120–350 mg/dl or for a mmol/L configured meter 117–351 mg/dl. For After low bG, the populated range is 050–100 mg/dl or for a mmol/L configured meter 54–100 mg/dl. For After meal, the populated range is 000–024 grams. <i>Dur Min – Dur Max</i> For After high bG, the range is 010000–060000. For After low bG, the range is 000500–003000. For After meal, the range is 010000–040000.

Roche Confidential Document
Released

Company Confidential

Document Number 7005079	Rev E
-----------------------------------	-----------------



Roche	CPDG A/C COMBO EXPERT DM
© 2009 All rights reserved	

8.6 Read Setup – [53] or ‘S’

Example: bG Reminder After high bG Enabled, Threshold = 200 mg/dl, bG Reminder Time Duration = 2 hours and 45 minutes.

```
→S
←S
→<TAB>
←<TAB>
→[74]
←[74]
→<TAB>
←<TAB>
→0
←0
→<CR>
←<ACK>
←<STX>{LEN}<TAB>1<TAB>200<TAB>024500<TAB>{CRC}<EOT>
→<ACK>
←<ACK>
```

Note 1: The bG Reminder Time Duration, for the After High bG and After Meal reminders, the minutes shall be a multiple of a quarter hour and the seconds value shall be zero. The bG Reminder Time Duration, for the After Low bG reminder, the minutes shall be a multiple of 5 minutes and the seconds value shall be zero.

Note 2: For the bG Reminder Threshold, the Meter shall return <TAB><TAB> if its value has not been set up.

Note 3: For the bG Reminder Time Duration, for the after High bG reminder, the Meter shall return <TAB><TAB> if its value has not been set up.

Roche Confidential Document
Released

Company
Confidential

Document Number	Rev
7005079	E



Roche © 2009 All rights reserved	CPDG A/C COMBO EXPERT DM
--	---------------------------------

8.6 Read Setup – [53] or ‘S’

S<TAB>[75] — Read Alarm Clock	The S<TAB>[75] command shall read the Selection and Time for a selected Alarm Clock.		
	Input Name Alarm Clock Number	Input Type Number	Input Range 0 – 7 0 = Alarm Clk 1 . . . 7 = Alarm Clk 8
	Output Name Selection	Output Type Number	Output Range 0 – 3 0 = Off 1 = bG Test 2 = Injection (Expert meters only) 3 = Other
	Time	Time	000000 – 234500
Example: Alarm Clock 3 set for bG Test at 8:15 AM →S ←S →<TAB> ←<TAB> →[75] ←[75] →<TAB> ←<TAB> →2 ←2 →<CR> ←<ACK> ←<STX>{LEN}<TAB>1<TAB>081500<TAB>{CRC}<EOT> →<ACK> ←<ACK> Note: For the Time, the minutes shall be a multiple of a quarter hour and the seconds value shall be zero. “000000” shall be returned for the Time if the time has not been set.			

**Company
Confidential**

Document Number	Rev
7005079	E



Roche © 2009 All rights reserved	CPDG A/C COMBO EXPERT DM
--	---------------------------------

S<TAB>[76] — Read Date Reminder	The S<TAB>[76] command shall return the Reminder Status, Time Set Status, Date, and Time for a selected Reminder Type.		
	Input Name Reminder Type	Input Type Number	Input Range 0 – 1 0 = Dr. Visit 1 = Lab Test
	Output Name Reminder Status Time Set Status Date Time	Output Type Number Number Date Time	Output Range 0 – 1 (0 = disabled) 0 – 1 0 = No 1 = Yes 030101 – 311231 000000 – 234500
	Example: Dr. Visit Reminder set for 3:45 PM, June 17, 2007 →S ←S →<TAB> ←<TAB> →[76] ←[76] →<TAB> ←<TAB> →0 ←0 →<CR> ←<ACK> ←<STX>{LEN}<TAB>1<TAB>1<TAB>070617<TAB>154500<TAB>{CR C}<EOT> →<ACK> ←<ACK> Note: For the Time, the minutes shall be a multiple of a quarter hour and the seconds value shall be zero. "000000" shall be return for the Time if the time has not been set.		

Roche Confidential Document
Released

Company Confidential

Document Number 7005079	Rev E
-----------------------------------	-----------------



Roche © 2009 All rights reserved	CPDG A/C COMBO EXPERT DM
--	---------------------------------

8.6 Read Setup – [53] or ‘S’

S<TAB>[77] — Read Infusion Set Reminder	The S<TAB>[77] command shall read the Status and Interval for the Infusion Set Reminder.		
	Output Name Status	Output Type Number	Output Range 0 – 1 (0 = disabled)
	Interval	Number	0 – 2 0 = 1 day 1 = 2 days 2 = 3 days
Example: Infusion Reminder Enabled and set for 2 days →S ←S →<TAB> ←<TAB> →[77] ←[77] →<CR> ←<ACK> ←<STX>{LEN}<TAB>1<TAB>1<TAB>{CRC}<EOT> →<ACK> ←<ACK> Note: This command is only intended for ACCU-CHEK® Combo meters.			

Roche Confidential Document
Released

Company Confidential

Document Number 7005079	Rev E
-----------------------------------	-----------------



Roche © 2009 All rights reserved	CPDG A/C COMBO EXPERT DM
--	---------------------------------

S<TAB>[78] — Read Bolus Advice Timeblock	The S<TAB>[78] command shall return the Time, Minimum and Maximum bG Values, Carb Ratio Insulin Value and Carb Ratio Carb Values, and Insulin Sensitivity Insulin Value and Insulin Sensitivity bG Values for a selected Bolus Advice Timeblock.		
	Input Name Bolus Advice Timeblock	Input Type Number	Input Range 0 – 7 0 = Timeblock 1 ... 7 = Timeblock 8
	Output Name End Time Minimum bG Value Maximum bG Value Carb Ratio Insulin Value Carb Ratio Carb Value Insulin Sensitivity Insulin Value Insulin Sensitivity bG Value	Output Type Time Number (mg/dl) Number (mg/dl) Float (units) Number (grams) Float (units) Number (grams)	Output Range 000000 – 233000 050 – 140 100 – 300 0.1 – 50.0 001 – 240 00.1 – 50.0 001 – 999

Roche Confidential Document
Released

Company Confidential

Document Number 7005079	Rev E
-----------------------------------	-----------------



Roche	CPDG A/C COMBO EXPERT DM
© 2009 All rights reserved	

8.6 Read Setup – [53] or ‘S’

**Example: Bolus Advice Timeblock 2: End Time = 4:30 PM,
Min bG = 30 mg/dl, Max bG = 250 mg/dl,
Carb Ratio Insulin Value = 4.5 units,
Carb Ratio Carb Value = 110 grams,
Insulin Sensitivity Insulin Value = 23.0 units,
Insulin Sensitivity bG Value = 800 mg/dl**

→S

←S

→<TAB>

←<TAB>

→[78]

←[78]

→<TAB>

←<TAB>

→1

←1

→<CR>

←<ACK>

←<STX>{LEN}<TAB>163000<TAB>030<TAB>250<TAB>04.5<TAB>11
0<TAB>23.0<TAB>800<TAB>{CRC}<EOT>

→<ACK>

←<ACK>

Note 1: If the Bolus Advice Timeblock Number corresponds to a timeblock that is not set up, the meter will NAK the command and update the Status Register with Invalid Parameter (F8H). The seconds value of the End Time is always set to zero.

Note 2: For the Carb Ratio Carb Value, the Meter shall return <TAB><TAB> if its value has not been set up.

Note 3: For the Insulin Sensitivity bG Value, the Meter shall return <TAB><TAB> if its value has not been set up.

Note 4: The Output Ranges for mmol/L configured meters are: Minimum bG Value 54-144 mg/dl and Maximum bG Value 100 – 270 mg/dl.

Roche Confidential Document
Released

Company
Confidential

Document Number	Rev
7005079	E



Roche © 2009 All rights reserved	CPDG A/C COMBO EXPERT DM
--	---------------------------------

8.6 Read Setup – [53] or ‘S’

S<TAB>[79] — Read Health Event Value	The S<TAB>[79] command shall read the Health Event Value for a selected Health Event.		
	Input Name Health Event Number	Input Type Number	Input Range 0 – 4 0 = Exercise 1 1 = Exercise 2 2 = Stress 3 = Illness 4 = Premenstrual
	Output Name Health Event Value	Output Type Signed Number (%)	Output Range “-50” – “+50” or “00”
	Example: Health Event Number 2, Stress = -13% →S ←S →<TAB> ←<TAB> →[79] ←[79] →<TAB> ←<TAB> →2 ←2 →<CR> ←<ACK> ←<STX>{LEN}<TAB>-13<TAB>{CRC}<EOT> →<ACK> ←<ACK>		

Roche Confidential Document
Released

Company Confidential

Document Number 7005079	Rev E
-----------------------------------	-----------------



Roche © 2009 All rights reserved	CPDG A/C COMBO EXPERT DM
--	---------------------------------

S<TAB>[80] — Read Bolus Advice Options	The S<TAB>[80] command shall return the Bolus Advice Options Meal Excursion (Meal Rise), Active Timeout, Offset Timeout, and Snack Limit values.		
	Output Name Meal Excursion (Meal Rise) Active Timeout Offset Timeout Snack Limit	Output Type Number (mg/dl) Time Time Number (grams)	Output Range 050 – 200 013000 – 080000 004500 – Active Timeout Value 00 – 24
	Example: Bolus Advice Meal Excursion = 70 mg/dl, Active Timeout = 3 hours 45 minutes, Offset Timeout = 2 hours 00 minutes, Snack Limit = 14 grams →S ←S →<TAB> ←<TAB> →[80] ←[80] →<CR> ←<ACK> ←<STX>{LEN}<TAB>070<TAB>034500<TAB>020000<TAB>14<TAB>{CRC}<EOT> →<ACK> ←<ACK> Note 1: For Active Timeout and Offset Timeout, the minutes shall be a multiple of a quarter hour and the seconds value shall be zero. Note 2: For the Snack Limit, the Meter shall return <TAB><TAB> if its value has not been set up.		

Roche Confidential Document
Released

Company Confidential

Document Number 7005079	Rev E
-----------------------------------	-----------------



Roche	CPDG A/C COMBO EXPERT DM
© 2009 All rights reserved	

8.6 Read Setup – [53] or ‘S’

S<TAB>[82] — Read Patient Hyper bG threshold	The S<TAB>[82] command shall read the Patient Hyper bG threshold.		
	Output Name Hyper bG Threshold	Output Type Number (mg/dl)	Output Range 180 – 350
	Example: Hyper bG Threshold = 200 mg/dl →S ←S →<TAB> ←<TAB> →[82] ←[82] →<CR> ←<ACK> ←<STX>{LEN}<TAB>200<TAB>{CRC}<EOT> →<ACK> ←<ACK> Note: The Output Range for a mmol/L configured meter is: 180- 351 mg/dl.		
S<TAB>[83] — Read Energy Units	The S<TAB>[83] command shall return the Energy Units setting.		
	Output Name Energy Units	Output Type Number	Output Range 0 – 2 0 = Cal 1 = KCal 2 = KJ
	Example: Energy Units = Calories →S ←S →<TAB> ←<TAB> →[83] ←[83] →<CR> ←<ACK> ←<STX>{LEN}<TAB>0<TAB>{CRC}<EOT> →<ACK> ←<ACK>		

Company
Confidential

Document Number	Rev
7005079	E



Roche © 2009 All rights reserved	CPDG A/C COMBO EXPERT DM
--	---------------------------------

8.6 Read Setup – [53] or ‘S’

S<TAB>[84] — Read Max Bolus Value	The S<TAB>[84] command shall read the Max Bolus Value.		
	Input Name Selection	Input Type Number	Input Range 0
	Output Name Max Bolus Value	Output Type Float (Insulin Units)	Output Range 00.0 – 50.0
	<p>Example: ICT Max Bolus Value = 23.5 units</p> <p>→S ←S →<TAB> ←<TAB> →[84] ←[84] →<TAB> ←<TAB> →0 ←0 →<CR> ←<ACK> ←<STX>{LEN}<TAB>23.5<TAB>{CRC}<EOT> →<ACK> ←<ACK></p> <p>Note: This command is only available on ACCU-CHEK® Expert meters. ACCU-CHEK® Combo meters do not support this command.</p> <p>Note 2: For the Max Bolus Value, the Meter shall return <TAB><TAB> if its value has not been set up.</p>		

Roche Confidential Document

Released

**Company
Confidential**

Page 36 of 49

Document Number
7005079

Rev
E



Roche © 2009 All rights reserved	CPDG A/C COMBO EXPERT DM
--	---------------------------------

8.7 Change Setup – [0C]			
Command Character [0C]	The [0C] commands shall allow the host to alter the settings of the meter.		
[0C]<TAB>1 — Set date	The [0C]<TAB>1 command shall set the meter's date.		
	Input Name Date	Input Type Date (YYMMDD)	Input Range 090101 – 311231
	Example: Set date to February 10, 2009 →[0C] ←[0C] →<TAB> ←<TAB> →1 ←1 →<CR> ←<ACK> →<STX>{LEN}<TAB>090210<TAB>{CRC}<EOT> ←<ACK> Note: Setting the meter to an invalid date will result in a meter error.		
[0C]<TAB>2 — Set time	The [0C]<TAB>2 command shall set the meter's time.		
	Input Name Time	Input Type Time (HHMMSS)	Input Range 000000 – 235959
	Example: Set time to 3:45 PM →[0C] ←[0C] →<TAB> ←<TAB> →2 ←2 →<CR> ←<ACK> →<STX>{LEN}<TAB>154500<TAB>{CRC}<EOT> ←<ACK> ←<ACK>		

Company Confidential

Document Number 7005079	Rev E
-----------------------------------	-----------------

Roche Confidential Document
Released

**8.7 Change Setup – [0C]****[0C]<TAB>7 — Set Language**

The [0C]<TAB>7 command shall set the selected language.

Input Name

Selected Language

Input Type

Number

Input Range

1 – 7, A, C, D, F, G

H, I, J, M, N, O, Q

1 = German

2 = English

3 = French

4 = Spanish

5 = Italian

6 = Dutch

7 = Swedish

A = Portuguese

C = Norwegian

D = Finnish

F = Danish

G = Czech

H = Hungarian

I = Slovenian

J = Russian

M = Polish

N = Slovakian

O = Brazilian

Portuguese

Q = Romanian

Example: Set Language to English (a stored language)

→[0C]

←[0C]

→<TAB>

←<TAB>

→7

←7

→<CR>

←<ACK>

→<STX>{LEN}<TAB>2<TAB>{CRC}<EOT>

←<ACK>

←<ACK>

Note: If the Selected Language Parameter is set to a number that does not correspond to any of the languages stored in the meter, the meter shall NAK the command and updates the Status Register with INVALID_PARAMETER (F8H).

**Company
Confidential**

Roche Confidential Document
Released



Roche © 2009 All rights reserved	CPDG A/C COMBO EXPERT DM
--	---------------------------------

8.7 Change Setup – [0C]

[0C]<TAB>9 — Set time format	The [0C]<TAB>9 command shall set the meter's time format.		
	Input Name Time Format	Input Type Number	Input Range 1 – 2 1 = 24 2 = 12
	Example: Set time format to 12-Hour →[0C] ←[0C] →<TAB> ←<TAB> →9 ←9 →<CR> ←<ACK> →<STX>{LEN}<TAB>2<TAB>{CRC}<EOT> ←<ACK> ←<ACK>		
[0C]<TAB>[71] — Set Key Click Status	The [0C]<TAB>[71] command shall set the key click status		
	Input Name Key Click Status	Input Type Number	Input Range 0 – 1 (0 = disable)
	Example: Disable Key Click →[0C] ←[0C] →<TAB> ←<TAB> →[71] ←[71] →<CR> ←<ACK> →<STX>{LEN}<TAB>0<TAB>{CRC}<EOT> ←<ACK> ←<ACK>		

Roche Confidential Document
Released

Company Confidential

Document Number 7005079	Rev E
-----------------------------------	-----------------



Roche	CPDG A/C COMBO EXPERT DM
© 2009 All rights reserved	

8.7 Change Setup – [0C]

[0C]<TAB>[72] — Set Beeper Volume	The [0C]<TAB>[72] command shall set the beeper volume.		
	Input Name Beeper Volume	Input Type Number	Input Range 0 – 3 0 = Off 1 = Low 2 = Medium 3 = High
	<p>Example: Set Beeper Volume to Low</p> <p>→[0C] ←[0C] →<TAB> ←<TAB> →[72] ←[72] →<CR> ←<ACK> →<STX>{LEN}<TAB>1<TAB>{CRC}<EOT> ←<ACK> ←<ACK></p> <p>Note: Due to safety concerns, the meter does not allow the user interface to disable both the vibrator and beeper. It is highly recommended that host applications ensure that the vibrator and beeper are not both disabled. The meter does not enforce this through the serial interface.</p>		

Roche Confidential Document
Released

Company Confidential

Document Number 7005079	Rev E
-----------------------------------	-----------------



Roche	CPDG A/C COMBO EXPERT DM
© 2009 All rights reserved	

8.7 Change Setup – [0C]			
[0C]<TAB>[73] — Set Vibrate Status	The [0C]<TAB>[73] command shall set the vibrate status.		
	Input Name Vibrate Status	Input Type Number	Input Range 0 – 1 (0 = disable)
	<p>Example: Set Vibrate status to Disabled.</p> <p>→[0C] ←[0C] →<TAB> ←<TAB> →[73] ←[73] →<CR> ←<ACK> →<STX>{LEN}<TAB>0<TAB>{CRC}<EOT> ←<ACK> ←<ACK></p> <p>Note: Due to safety concerns, the meter does not allow the user interface to disable both the vibrator and beeper. It is highly recommended that host applications ensure that the vibrator and beeper are not both disabled. The meter does not enforce this through the serial interface.</p>		

Roche Confidential Document
Released

Company Confidential

Document Number 7005079	Rev E
-----------------------------------	-----------------



Roche © 2009 All rights reserved	CPDG A/C COMBO EXPERT DM
--	---------------------------------

8.8 Enable/Disable Timers – [5A] or ‘Z’			
Command Character [5A] or ‘Z’	The 5A or Z command shall toggle the meter's timeout timers.		
Z<TAB>0 — Disable selected timeouts	The Z<TAB>0 command shall disable the selected timeouts.		
	Input Name Timers	Input Type Function	Input Range ‘0’ – ‘9’, ‘A’ – ‘F’
	<p>Example: Disable all timeouts</p> <p>→Z ←Z →<TAB> ←<TAB> →0 ←0 →<TAB> ←<TAB> →F ←F →<CR> ←<ACK> ←<ACK></p> <p>Note: The ‘0’ in this command is the digit zero and not the letter ‘O’.</p> <p>The input parameter is a bit mapped of the four available timeouts. The bits are mapped as follows:</p> <p>0 - Disable auto-power off timer 1 - Disable command timeout 2 - Disable data exchange timeout 3 - Disable inter-character timeout</p>		

Roche Confidential Document
Released

Company Confidential

Document Number 7005079	Rev E
-----------------------------------	-----------------



Roche	CPDG A/C COMBO EXPERT DM
© 2009 All rights reserved	

8.8 Enable/Disable Timers – [5A] or ‘Z’

Z<TAB>1 — Enable selected timeouts	The Z<TAB>1 command shall enable selected timeouts.		
	Input Name Timers	Input Type Function	Input Range ‘0’ – ‘9’, ‘A’ – ‘F’
	<p>Example: Enable all timeout</p> <p>→Z ←Z →<TAB> ←<TAB> →1 ←1 →<TAB> ←<TAB> →F ←F →<CR> ←<ACK> ←<ACK></p> <p>Note 1: The input parameter is a bit mapped of the four available timeouts. The bits are mapped as follows:</p> <p>Bit 0 - Enable auto-power off timer Bit 1 - Enable command timeout Bit 2 - Enable data exchange timeout Bit 3 - Enable inter-character timeout</p> <p>Note 2: Auto-power off timeout – If the meter goes more than 120 seconds without a button, strip or command being sent it automatically shuts down.</p>		

Roche Confidential Document
Released

Company Confidential

Document Number 7005079	Rev E
-----------------------------------	-----------------



Roche © 2009 All rights reserved	CPDG A/C COMBO EXPERT DM
--	---------------------------------

8.9 Obtain Number of Results – [60]

Command Character [60]	The [60] command shall return the number of results stored in the meter.		
	Output Name Number Records	Output Type Number	Output Range 0 – 1100
	Example: Number of results = 1100 →[60] ←[60] →<CR> ←<ACK> ←<STX>{LEN}<TAB>1100<TAB>{CRC}<EOT> →<ACK> ←<ACK> Note: Corrupted, Blank, and pump data only records are not transmitted, however are included in the total number of records.		

8.10 Send Results from Start to End – [61] or ‘a’

Command Character [61] or ‘a’	The [61] or ‘a’ command shall return the requested bG results, cG results, carb data, and insulin data.		
	Input Name Start Value End Value	Input Type Number Number	Input Range 1 – Number of records stored 1 – Number of records stored
	Output Name Glucose value Time Date Flags DM Data Block	Output Type Number (mg/dl) or empty Time (HHMM) Date (YYMMDD) Fixed Hex Double Word Empty or see definition below	Output Range 0 – 999 0000 – 2359 010101 – 311231 See bit pattern below. Empty or see definition below

**Company
Confidential**

Document Number	Rev
7005079	E

**8.10 Send Results from Start to End – [61] or ‘a’****Example: Read 10 results starting with result 1**

```
→[61]
←[61]
→<TAB>
←<TAB>
→1
←1
→<TAB>
←<TAB>
→1
←1
→0
←0
→<CR>
←<ACK>
←<STX>{LEN}<TAB>120<TAB>2359<TAB>030612<TAB>00000010
<TAB><TAB>{CRC}<ETX>
→<ACK>
...
←<STX>{LEN}<TAB>120<TAB>1234<TAB>030612<TAB>00000020
<TAB><TAB>{CRC}<EOT>
→<ACK>
←<ACK>
```

Note 1: The flags field will transmit all zeros if no flags are set.

Note 2: Corrupted, blank, or pump data only records are not transmitted, however are included in the total number of records.

Note 3: For bG records, if the Glucose Value does not exist, the [61] command shall return an empty Glucose Value (<TAB><TAB>).

Note 4: For cG records, the [61] command shall return an empty DM Data Block.

Note 5: bG and cG results shall be transmitted in the reverse order in which they were stored.

**Company
Confidential**



Roche

© 2009 All rights reserved

CPDG A/C COMBO EXPERT DM**8.10 Send Results from Start to End – [61] or ‘a’**

	DM Data Block:	
	Number of Data Blocks	If the record contains DM data blocks, the first field shall be the number of data blocks that will follow.
	Meal Time Information	If the record contains Meal Time information, the DM data block shall contain <TAB>37<TAB>event code<TAB><TAB><TAB><TAB> where <i>event code</i> is 1 for premeal, 2 for post meal, 4 for bedtime, 74 for other.
	Health Event Information	If the record contains Health Event information, the DM data block shall contain <TAB>37<TAB>health event<TAB><TAB><TAB><TAB> where health event is 75 for Exercise 1, 76 for Exercise 2, 29 for Stress, 73 for Premenstrual, 20 for Fasting, and 31 for Illness.
	Meal Time and Health Event Information	If both Meal time and Health Event information exists for a record, the DM data block shall contain <TAB>37<TAB>event code<TAB><health event><TAB><TAB><TAB> where event code is 1 for premeal, 2 for post meal, 4 for bedtime, 74 for other and where health event is 75 for Exercise 1, 76 for Exercise 2, 29 for Stress, 73 for Premenstrual, 20 for Fasting, and 31 for Illness.

**Company
Confidential**

Page 46 of 49

Document Number

7005079

Rev

E



Roche © 2009 All rights reserved	CPDG A/C COMBO EXPERT DM
-------------------------------------	--------------------------

8.10 Send Results from Start to End – [61] or ‘a’

	Carb Amount Information	If the record contains Carbohydrate information, the DM data block shall contain <TAB>12<TAB>value<TAB> where value is 3 characters representing the amount of carbohydrates in grams.
	Pen/Syringe Insulin Bolus	If the record contains Insulin Bolus information, the DM data block shall contain <TAB>34<TAB>1<TAB>29<TAB>insulin dosage<TAB> where insulin dosage is in the Float format (1 to 4 characters containing a decimal point).
	Pen/Syringe Insulin Basal	If the record contains Insulin Basal information, the DM data block shall contain <TAB>34<TAB>2<TAB>30<TAB>insulin dosage<TAB> where insulin dosage is in the Float format (1 to 4 characters containing a decimal point).
<p>Example of a result record with a DM data block:</p> <p>bG result was 120 mg/dl; it occurred at 11:59 PM on June 12, 2003, Bed time event; Health Event is Illness, Pen Syringe Insulin bolus with an amount of 4.5 IU.</p> <p>Result Packet: <STX>{LEN}<TAB>120<TAB>2359<TAB>030612<TAB>00000000<TAB>2<TAB>37<TAB>4<TAB>31<TAB><TAB><TAB>34<TAB>1<TAB>0<TAB>4.5<TAB>{CRC}<ETX></p>		
Flags:		
0x00000000	No flags	
0x00000004	Result too low	
0x00000008	Result too high	
0x00000010	Control solution level 2	

**Company
Confidential**

Page 47 of 49

Document Number
7005079

Rev
E

Roche Confidential Document
Released



Roche © 2009 All rights reserved	CPDG A/C COMBO EXPERT DM
--	---------------------------------

8.10 Send Results from Start to End – [61] or 'a'

	0x00000020	Control solution level 1
	0x00000040	Result out of temperature range
	0x00000200	Result below hypo
	0x00000800	User's result below his personal target range or control result below control's target range
	0x00001000	User's result above his personal target range or control result above control's target range
	0x00008000	Control not identified
	0x00020000	Result above hyper
	0x00040000	Before Meal (Flag not used in all models)
	0x00080000	After Meal (Flag not used in all models)

8.11 Reset Results Memory – [52] or 'R'

Command Character [52] or 'R'	The [52] or 'R' command shall clear results memory.
	Example: Clear results memory →R ←R →<CR> ←<ACK> ←<ACK> Note: This function clears all results memory. This includes bG results, cG results, and Bolus Advice history.

**Company
Confidential**

Page 48 of 49

Document Number	Rev
7005079	E



Roche © 2009 All rights reserved	CPDG A/C COMBO EXPERT DM
--	---------------------------------

9 Status Register Values

Status Value (hex)	Status Register Value Meaning
0000h	No Errors
0001h – 00EFh	Internal Meter Errors
00F0h	Command Canceled
00F1h	STX Expected Error
00F2h	Length Expected Error
00F3h	Not used.
00F4h	Not Used
00F5h	Not used.
00F6h	IR Data Overrun
00F7h	Invalid Number of Bytes
00F8h	Invalid Parameter
00F9h	Invalid Number of Parameters
00FAh	Receive Buffer Full
00FBh	Communication Timeout
00FCh	Command Not Implemented
00FDh	Command Aborted
00FEh	Not Valid Command
00FFh	Initial Communication
0100h – 0430h	Internal Meter Errors
0431h	Bolus Advice Setup Error
0432h – 0440h	Internal Meter Errors
0441h	Setup Wizard Setup Error
0442h – FFFFh	Internal Meter Errors

**Company
Confidential**

Page 49 of 49

Document Number
7005079

Rev
E